<국내대학원 석사야간과정 훈련과제 연구보고서>

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PROTOCOL FOR ANIMAL TISSUE

Before experiment

- Prepare the water bath to 56°C.
- Equilibrate buffer AE to room temperature.
- All centrifugation should be performed at room temperature.
- Buffer TB may precipitate at cold ambient temperature.
- If so, completely dissolve it in 37°C water bath.

I. Disrupt 20 mg of tissue as described on step Ia, Ib or Ic, depending on the sample type.

To disrupt the sample finer will accelerate lysis and reduce the lysis time.

For spleen tissue, up to 10 mg can be processed.

If the starting sample is larger than 25 mg (if spleen, 10 mg), increase the volume of buffer TL proportionally. For 50 mg of liver tissue, 400 ul of buffer TL is required.

For over 50 mg of tissue, Exgene™ Tissue SV Midi or MAXI is recommended.

Ia. For soft tissue, such as liver or brain, put 20 mg of the tissue into 1.5 ml tube, add 200 ul of Buffer TL, and homogenize thoroughly with microhomogenizer.

Ib. If microhomogenizer is not available or the tissue is not soft, grind the tissue to a fine powder with liquid nitrogen in a pre-chilled mortar and pestle. Put 20 mg of the powdered tissue into 1.5 ml tube. Add 200 ul of Buffer TL and pulse-vortex for 15 seconds.

Ic. If neither Ia nor Ib is available, mince 20 mg of tissue with sharp blade or scalpel as small as possible. Put the tissue into a 1.5 ml tube. Add 200 ul of Buffer TL and pulse-vortex for 15 seconds.

***** Alternatively, tissue samples can be effectively disrupted using some instruments, such as a rotor-stator homogenizer or a bead-beater.**

2. Add 20 ul of Proteinase K solution. Mix completely by vortexing or pipetting. Incubate at 56°C until the sample is completely lysed.

It is essential to mix the components completely for proper lysis.

If the sample amount is larger than 25 mg (if spleen, 10 mg), increase the amount of Proteinase K proportionally. For 50 mg of liver tissue, 40 ul of Proteinase K solution is required.

Lysis time varies from 10 min to 3 hr usually depending on the type of tissue processed and the disruption method (step1). The lysate should become translucent without any particles after complete lysis. Overnight lysis does not influence the preparation. If the sample tube is incubated in water bath or heating block, vortex occasionally (2-3 times per hour) during incubation to lyse readily. *Lysis in shaking water bath, shaking incubator or agitator would be best for efficient lysis.*

3. Check!! If Buffer TB precipitates, pre-heat in a 37°C water bath to dissolve completely.

Precipitated buffer will cause significant decrease in recover yield.

4. Spin down the tube briefly to remove any drops from inside of the lid.

5. Optional: If RNA-free DNA is required, add 4 ul of RNase solution (100 mg/ml, Cat.No.117-960), vortex to mix thoroughly, and incubate for 2 min at room temperature.

GeneAll® Exgene™ Tissue SV column has stronger affinity to DNA than RNA. Although the treatment of RNase is omitted, RNA occupies very small portion of eluates. RNA may inhibit some downstream enzymatic reactions, but will not inhibit PCR.

6. Add 400 ul of Buffer TB. Immediately vortex the tube to mix thoroughly. Spin down the tube briefly to remove any drops from inside of the lid.

If the sample is larger than 25 mg (if spleen, 10 mg), increase the volume of buffer TB proportionally. For 50 mg of liver tissue, 800 ul of buffer TB is required.

7. Apply the mixture to the SV column. Centrifuge for 1 min at 6,000 xg above (>8,000 rpm). Replace the collection tube with new one (provided).

If more than 25 mg (10 mg spleen) of tissue is processed, apply the mixture twice; apply 700 ul of the mixture, spin down, discard the filtrate, re-insert to the empty collection tube, and repeat the step again until all of the mixture is applied to the SV column. If the mixture has not passed completely through the membrane, centrifuge again at full speed until all of the solution has been passed through.

Centrifuge at maximum speed will not affect the DNA recovery.

8. Add 600 ul of Buffer BW. Centrifuge for 30 sec at 6,000 xg above (>8,000 rpm). Replace the collection tube with new one (provided).

9. Apply 700 ul of Buffer TW. Centrifuge for 30 sec at 6,000 xg above (>8,000 rpm). Discard the filtrate and reinsert the SV column back into the collection tube.

10. Centrifuge at full speed (above 13,000 xg) for 1 min to remove residual wash buffer. Place the SV column in a fresh 1.5 ml tube (not provided).

Care must be taken at this step for eliminating the carryover of buffer TW. If carryover of buffer TW occurs on column membrane, centrifuge again for 1 min at full speed before transferring to the new 1.5 ml tube.

2. PCR(polymerase chain reaction) & Electrophoresis

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
: Total 192 samples, 96 lung & 96 hilar LN

나. Preparation of a set of primer for PCV3 detection

TABLE 1 List of primers used in this study

Primer name	Nucleotide sequence	Primer location (nt) ^a	Product length (bp)	Purpose
PCV3-1-F	T T A C T T A G A G A A C G G A C T T G T A A C G	1339-1363	649	Detection
PCV3-1-R	A A A T G A G A C A C A G A G C T A T A T T C A G	1987-1965		

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3. Multiplex real-time PCR for pcv2 and pcv3

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Reference

1. Allan, G., Krakowka, S., Ellis, J., Charreyre, C., 2012. Discovery and evolving history of two genetically related but phenotypically different viruses, porcine circoviruses 1 and 2. *Virus Res.* 64, 4-9.
2. Hye-Ryung Kima, Yu-Ri Parka, Da-Rae Lima, Min-Ji Parka, Ji-Young Parkb, Seong-Hee Kimb, Kyoung-Ki Leeb, Young S. Lyooc, Choi-Kyu Park, Multiplex real-time polymerase chain reaction for the differential detection of porcine circovirus 2 and 3. *Journal of Virological Methods* 250 (2017) 11-16
3. Ku, X., Chen, F., Li, P., Wang, Y., Yu, X., Fan, S., Qian, P., Wu, M., He, Q., 2017. Identification and genetic characterization of porcine circovirus type 3 in China. *Transbound. Emerg. Dis.* 64, 703-708.
4. Palinski, R., Pi~neyro, P., Shang, P., Yuan, F., Guo, R., Fang, Y., . . . Hause, B. M. (2016). A novel porcine circovirus distantly related to known circoviruses is associated with porcine dermatitis and nephropathy syndrome and reproductive failure. *Journal of Virology*, 91(1), e01879-16.
5. Seung-Chai Kim¹, Salik Nazki¹, Sam Kwon¹, Jhung-Hyun Juhng¹, Kyeong-Hwan Mun¹, Do-Young Jeon¹, Chang-Gi Jeong¹, Amina Khatun¹, Seog-Jin Kang² and Won-Il Kim^{1*}, Kim et al. *BMC Veterinary Research* (2018) 14:294
6. Tischer, I., Gelderblom, H., Vettermann, W., Koch, M.A., 1982. A very small porcine virus with circular single-stranded DNA.

Nature 295, 64-66.

7. Kwon, T., Lee, D.U., Yoo, S.J., Je, S.H., Shin, J.Y., Lyoo, Y.S., 2017a. Genotypic diversity of porcine circovirus type 2 (PCV2) and genotype shift to PCV2d in Korean pig population. *Virus Res.* 228, 24-29.
8. Opriessnig, T., Meng, X.J., Halbur, P.G., 2007. Porcine circovirus type 2 associated disease: update on current terminology, clinical manifestations, pathogenesis, diagnosis, and intervention strategies. *J. Vet. Diagn. Invest.* 19, 591-615.

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Appendix 4. Comparison of Multiplex quantitative real-time PCR